

Chemiluminescence Immunoassay Solution

We have some Special Items!

- Fully automatic operation
- High-throughput sample processing
- Small size: flexible and convenient
- Accurate results
- High sensitivity



Beijing Diagreat Biotechnologies Co.,Ltd

Diagreat Automatic Chemiluminescence Immunoassay Analyzer Therapeutic Drug Monitoring (TDM) Solution

Diagreat Automatic Chemiluminescence Immunoassay DF200i and DF100i analyzer adopts chemiluminescent assay of alkaline phosphatase- labeled detection technology.

The compatible reagents are independently developed from raw materials to products. The monoclonal antibodies can specifically identify small molecule drugs. The small molecule directional conjugated antigens can specifically modify. And the reagents can identify small molecule directional coupling technology.

Covering 20 TDM items, Diagreat keeps pace with other IVD leaders. Diagreat spares no efforts in constantly increasing test items based on medical demands.

DF200i Automatic Chemiluminescence Immunoassay Analyzer (Desktop Type)



Specifications Output	LCD Screen Internal printer LIS/HIS Connection
Memory	8GB
Size	865 × 620 × 1250mm
Weight	129KGS
Power voltage	100V-240V~, 50Hz/60Hz
Advantages	Test speed: 200 Tests/hour. Test time for the first sample: ≤15min. Test tube locations: 25. Sample locations: 60, disk system support continuous loading. Calibration type: two-point or three-point calibration. Operation interface: English. Multi-sample types: serum, plasma and urine.

DF100i Automatic Chemiluminescence Immunoassay Analyzer (Portable Type)



Specifications Output	Touch screen External printer LIS Connection
Memory	8GB
Size	433 × 678 × 638mm
Weight	62KGS
Power voltage	100-240V~, 50/60Hz
Advantages	Test speed: 120 Tests/hour. Test time for the first sample: ≤15min. Test tube locations: 10. Sample locations: 30, disk system support continuous loading. Calibration type: two-point or three-point calibration. Operation interface: English. Multi-sample types: serum, plasma and urine

Reagents Kits List

TDM Chemicals

Special Items!

MPA	Mycophenolic Acid	TEIC	Teicoplanin
MTX	Methotrexate	tRSP	Total Risperidone
VM	Vancomycin	PTX	Paclitaxel
Carb	Carbamazepine	LTG	Lamotrigine
PHT	Phenytoin	Vorl	Voriconazole
VPA	Valproic Acid	5-FU	5-Fluorouracil
Theo	Theophylline	PMB	Polymyxin B
DIG	Digoxin	RIF	Rifampicin
CsA	Cyclosporine A	QTP	Quetiapine
FK506	Tacrolimus	PHB	Phenobarbital
Gent	Gentamicin	DOC	Docetaxel
LZD	Linezolid	IMA	Imatinib
MEM	Meropenem	Tarl	Total aripiprazole
RAPA	Rapamycin		

Vitamins

25-OH-VD	Total 25-hydroxy Vitamin D
VB12	Vitamin B12
FA	Folic Acid

Cardiac

D-Dimer	D-Dimer
NT-proBNP	N-terminal pro-B-type Natriuretic Peptide
S100-β	S100-β Protein
hs-cTnI	High-sensitivity Troponin I
hs-cTnT	High-sensitivity Troponin T
BNP	B-type Natriuretic Peptide
MYO	Myoglobin
CK-MB	Creatine Kinase Isoenzymes MB
cTnI	Troponin I
cTnT	Troponin T
HFABP	H-Fatty acid binding protein
LP-PLA2	Lipoprotein Phospholipase A2

Anemia

Fer	Ferritin
FA	Folic Acid

Immunological

IgE	Immunoglobulin E
-----	------------------

TDM Biologics

Special Items!

Adalimumab	Cetuximab
Free Anti-adalimumab	Anti-Cetuximab
Infliximab	Neutralizing Anti-Cetuximab
Free Anti-infliximab	Trastuzumab
Golimumab	Anti-Trastuzumab
Free Anti-golimumab	Neutralizing Anti-Trastuzumab
Cetuximab	Ustekinumab
Free Anti-cetuximab	Free Anti-Ustekinumab
Bevacizumab	Neutralizing Anti-Ustekinumab
Free Anti-bevacizumab	Vedolizumab
Rituximab	Free Anti-Vedolizumab
Free Anti-rituximab	Neutralizing Anti-Vedolizumab
Obinutuzumab	Bevacizumab
Free Anti-Obinutuzumab	Free Anti-Bevacizumab
Adalimumab	Neutralizing Anti-Bevacizumab
Anti-Adalimumab	Rituximab
Neutralizing Anti-Adalimumab	Free Anti-Rituximab
Infliximab	Neutralizing Anti-Rituximab
Anti-Infliximab	Obinutuzumab
Total Anti-Infliximab	Free Anti-Obinutuzumab
Neutralizing Anti-Infliximab	Neutralizing Anti-Obinutuzumab
Golimumab	Etanercept
Anti-Golimumab	Free Anti-Etanercept
Neutralizing Anti-Golimumab	Neutralizing Anti-Etanercept

Neonatal

TSH	Thyroid stimulating hormone / Thyrotropin
17OHP	17α-hydroxyprogesterone
CK-MM	Neonatal Creatine Kinase-MM

Hormones

E2	Estradiol
FSH	Follicle Stimulating Hormone
Prog	Progesterone
TESTO	Testosterone
β-HCG	β-Human Chorionic Gonadotropin
PRL	Prolactin
LH	Luteinizing Hormone

Tumor Marker

CA19-9	Carbohydrate Antigen 19-9
CA72-4	Cancer Antigen 72-4
CA50	Cancer Antigen 50
CA125	Cancer Antigen 125
CA15-3	Cancer Antigen 15-3
CA242	Cancer Antigen 242
CEA	Carcinoembryonic Antigen
AFP	Alpha-fetoprotein
HE4	Human Epididymis Protein 4
CYFRA21-1	Cytokeratin 19 Fragment 21-1
NSE	Neuron Specific Enolase
ProGRP	Gastrin Releasing Peptide Precursor
tPSA	Total prostate specific antigen
fPSA	Free prostate specific antigen
SCCA	Squamous Cell Carcinoma Antigen

Thyroid

T3	Total Triiodothyronine
T4	Total Thyroxine
PTH	Parathyroid Hormone
TSH	Thyrotropin
FT3	Free Triiodothyronine
FT4	Free Thyroxine
Tg	Thyroglobulin
Anti TPO	Anti-thyroid Peroxidase Antibody
Anti-Tg	Anti-thyroglobulin

Hepatobiliary

LN	Laminin
CIV	Collagen IV
PIIINP	N-Terminal Procollagen III Propeptide
HA	Hyaluronic Acid

Inflammation

PCT	Procalcitonin
SAA	Serum Amyloid A
IL-2R	Interleukin-2 Receptor
IL-6	Interleukin-6
IL-8	Interleukin-8
IL-1 β	Interleukin-1 β
IL-10	Interleukin-10
CRP	C Reactive Protein
TNF- α	α Tumor Necrosis Factor

Gastric

PG I	Pepsinogen I
PG II	Pepsinogen II
G17	Gastrin 17

Glucose metabolism

CP	C-Peptide
INS	Insulin

Infectious Disease

HBsAg
Anti-HBs
HBeAg
Anti-HBe
Anti-HBc
HCV Ab/Ag Combi
HIV Ab/Ag Combi
Syphilis Ab (TP)

To learn more about the DIAGREAT CLIA KITS please contact your local representative!

For in vitro diagnostics use unless otherwise noted. Not all products are available in all territories due to different national regulations.

E-mail: info@diagreat.com
export@diagreat.com

Whatsapp/Mob:
+86-13524602677
+86-18500563281

[Http://www.diagreat.com](http://www.diagreat.com)



Add: Floor 1-3, Building 19, No. 2 Yard, Huanke
Middle Road, Tongzhou District, Beijing, China.